


ASSUFIL®



ARTI KODE PADA KEMASAN:



COMPOSITION	Polymer of glycolic acid
DESCRIPTION	Braided and coated synthetic absorbable suture (monofilament for size USP 100 - 8/0) EP 5/2 - 5/8
COLOR	White (SAC Mark 2 - C1, A0720) or uncolored
SIZE AND STANDARD NEEDLES	EP from 0,2 to 5, USP from 100 to 5. Needle is attached to 300 (30) or 450 (45) stainless steel absorbable needles. "knot-tethered" needle type is also available to allow the detachment of the needle by a simple movement of the needle holder.
TENSILE STRENGTH	After 7 days from implant - 90% of the tensile strength remains, after 14 days - 75%, after 21 days - 50% and after 28 days - 30%.
Absorption	Complete within 60-90 days.
INDICATIONS	In general surgery, thoracic surgery, orthopedics, gynecology, reproductive surgery, neurosurgery, in soft tissue approximation, when it is necessary an absorbable surgical suture.

THREAD LENGTH	SUFFIX
cm 45	AH
cm 60	AM
cm 70	AP
cm 75	AQ
cm 90	AT
cm 100	AV



Exclusive Distributor:


PT. LIMA MITRA ANDALAN
 Ruko Grand Caman
 Jl. Caman Raya No.54, Jati Bening, Bekasi
 Indonesia 17412
 Phone : +62 21 849 777 40 Fax : +62 21 849 0 1984
 E-mail : info@limamitra.com web : www.limamitra.com

ASSUFIL®

ABSORBABLE BRAIDED

BENANG

PGA - Polyglycolic Acid



Efektifitas

- ⊙ Dapat memegang jaringan selama masa penyembuhan (masa pemulihan jaringan 42 hari)
- ⊙ Simpul kuat, sehingga jaringan tidak mudah lepas

Mudah digunakan

- ⊙ Lembut, karena terbuat dari struktur benang yang halus dengan pelapis yang tipis (kurang dari 0.5%), sehingga mudah saat menarik benang pada waktu penjahitan
- ⊙ Dengan teknologi laser, penyatuan benang dan jarum menjadi lebih sempurna, sehingga kuat dan tidak menimbulkan trauma jaringan
- ⊙ Akses langsung, hanya dengan 2 langkah

Nyaman bagi pasien

- ⊙ Hasil penyembuhan luka operasi sangat baik
- ⊙ Benang tidak perlu diangkat, sehingga tidak sakit dan tidak membuat cemas pasien

Aman

- ⊙ Respon inflamasi minimal, karena penyerapan secara hidrolisis
- ⊙ Masa absorpsi hari ke 60-90

JARUM

Berkualitas tinggi













- ⊙ Terbuat dari stainless steel dengan kualitas AISI 400. Kualitas stainless steel setara dengan jarum benang bedah cardiovascular
- ⊙ Jarum tidak mudah bengkok

Ketajaman

- ⊙ Tetap tajam walaupun dilakukan sampai dengan 10x penetrasi
- ⊙ Dilapisi silikon

Indikasi

- ⊙ Bedah Umum
- ⊙ Bedah Digestive
- ⊙ Bedah Gynecology
- ⊙ Bedah Urology
- ⊙ Bedah Orthopedic

AGO -NEEDLE-AGUJA mm.	INDICATION	COLOR	THREAD X CM	USP	CODE	PACKS	UNIT PRICE (PCS)
 TAPER 1/2 CIRCLE							
 HR 17,4	Urologi, Usus	Violet	70	4 - 0	FW375	36	Contact us
 HR 21,8	Urologi, Usus	Violet	70	3 - 0	FV389	36	Contact us
 HR 25,9	Otot	Violet	70	3 - 0	FV386	36	Contact us
	Otot, Subcutis	Violet	70	2 - 0	FU395	36	Contact us
 HRG 36,4	Subcutis	Violet	70	2 - 0	FU440	36	Contact us
 HRG 39,9	Otot, Fascia	Violet	90	1	FR458	36	Contact us
 HRG 48	Uterus	Violet	90	1	FR479	36	Contact us
	Uterus	Violet	90	2	FQ478	36	Contact us
 CUTTING 3/8 CIRCLE							
 TS 13,1	Kulit	Violet	45	6 - 0	FY485	36	Contact us
	Kulit	Violet	45	5 - 0	FX503	36	Contact us
 TS 16,2	Kulit	Violet	45	5 - 0	FX499	36	Contact us
	Subcutis	Violet	45	4 - 0	FW500	36	Contact us
 TS 18,7	Subcutis	Violet	70	4 - 0	FW515	36	Contact us
	Subcutis	Violet	70	3 - 0	FV516	36	Contact us
 TS 24,3	Subcutis	Violet	70	3 - 0	FV533	36	Contact us
	Subcutis	Violet	70	2 - 0	FU534	36	Contact us

Harga berlaku mulai bulan Januari 2016

Harga belum termasuk PPN

Harga dapat berubah sewaktu - waktu tanpa pemberitahuan terlebih dahulu