


ASSUPRO®



ARTI KODE PADA KEMASAN:



COMPOSITION	Polypropylene.
DESCRIPTION	Synthetic non-absorbable monofilament suture.
COLOR	Blue (C.I. 74160 phthalocyanine).
SIZES AND STANDARD NEEDLES	EP from da 0,2 to 4. USP from 10/0 to 1. Assupro is attached to AISI 300 or AISI 400 stainless steel atraumatic needles. It is also packed in pouches containing teflon pledget or "sundry accessories" (polypropylene beads, collars and buttons) used to fix the thread end in intradermal sutures.
INDICATIONS	In cardiovascular surgery, neurosurgery, general surgery, ophthalmic surgery, plastic surgery, skin surgery, thoracic surgery, when it is necessary to suture slow healing tissues or tissues requiring a non absorbable suture. As most monofilament sutures, Assupro is indicated for suturing infected or contaminated tissues.

THREAD LENGTH	SUFFIX
cm 45	AH
cm 60	AM
cm 70	AP
cm 75	AC
cm 90	AT
cm 100	AV



ASSUPRO®

NON - ABSORBABLE MONOFILAMENT



BENANG

Polipropilene

Efektifitas

- ⊙ Dapat memegang jaringan selama masa penyembuhan (masa pemulihan jaringan 14 hari)
- ⊙ Simpul kuat, sehingga jaringan tidak mudah lepas
- Mudah digunakan**
- ⊙ Lembut, karena terbuat dari struktur benang yang halus dengan pelapis yang tipis (kurang dari 0.5%), sehingga mudah saat menarik benang pada waktu penjahitan
- ⊙ Dengan teknologi laser, penyatuan benang dan jarum menjadi lebih sempurna, sehingga kuat dan tidak menimbulkan trauma jaringan
- ⊙ Akses langsung, hanya dengan 2 langkah
- Nyaman bagi pasien**
- ⊙ Hasil penyembuhan luka operasi sangat baik

JARUM

Berkualitas tinggi

- ⊙ Terbuat dari stainless steel dengan kualitas AISI 400. Kualitas stainless steel setara dengan jarum benang bedah cardiovascular Jarum tidak mudah bengkok












Ketajaman









- ⊙ Tetap tajam walaupun dilakukan sampai dengan 10x penetrasi
- ⊙ Dilapisi silikon

Indikasi

- ⊙ Bedah Umum
- ⊙ Bedah Digestive
- ⊙ Bedah Gynecology
- ⊙ Bedah Urology
- ⊙ Bedah Orthopedic
- ⊙ Bedah Thorax
- ⊙ Bedah Syaraf
- ⊙ Bedah Jantung
- ⊙ Bedah Mata
- ⊙ Bedah Plastic

Exclusive Distributor:

AGO - NEEDLE-AGUJA mm.	INDICATION	COLOR	THREAD X CM	USP	CODE	PACKS	UNIT PRICE (PCS)
 TAPER 1/2 CIRCLE FLAT BODY							
	Usus	Blue	75	4 - 0	PW379	36	Contact us
			75	3 - 0	PV378	36	Contact us
			75	2 - 0	PU395	36	Contact us
	Hernia	Blue	75	0	PS436	36	Contact us
			75	1	PR458	36	Contact us
 TAPER 3/8 CIRCLE FLAT BODY							
	Cardiovascular	Blue	75	8 - 0	PA875K	36	Contact us
			75	7 - 0	PZ925K	36	Contact us
	Cardiovascular	Blue	75	7 - 0	PZ890	36	Contact us
			75	6 - 0	PY891	36	Contact us
	vascular	Blue	75	7 - 0	PZ892K	36	Contact us
			75	6 - 0	PY893K	36	Contact us
			75	5 - 0	PX894K	36	Contact us
 TAPER 1/2 CIRCLE ULTRA							
	Cardiovascular	Blue	60	6 - 0	PY87AKAT	36	Contact us
	Cardiovascular	Blue	90	5 - 0	PX89AK	36	Contact us
			90	4 - 0	PW90AK	36	Contact us
	Cardiovascular	Blue	90	4 - 0	PW92AK	36	Contact us
			90	3 - 0	PV93AK	36	Contact us
			90	2 - 0	PU94AK	36	Contact us

AGO - NEEDLE-AGUJA mm.	INDICATION	COLOR	THREAD X CM	USP	CODE	PACKS	UNIT PRICE (PCS)
 CUTTING 3/8 CIRCLE							
	Kulit	Blue	75	6 - 0	PY498	36	Contact us
			75	5 - 0	PX499	36	Contact us
	Kulit	Blue	75	4 - 0	PW515	36	Contact us
			75	3 - 0	PV516	36	Contact us
	Kulit	Blue	75	3 - 0	PV533	36	Contact us
	Kulit	Blue	75	2 - 0	PU543	36	Contact us
	Kulit	Blue	75	2 - 0	PU539	36	Contact us
 TAPER 3/8 CIRCLE FLAT BODY							
	Cardiovascular	Blue	2x75	7 - 0	PZ890S	36	Contact us