


GLICOFIL® LAC



ARTI KODE PADA KEMASAN:



COMPOSITION	PGLA Co-polymer of glycolide (90%) and lactide (10%)
DESCRIPTION	Braided and coated synthetic absorbable suture (monofilament for size: USP 10/0 - 9/0; EP 0.2 - 0.3)
COLOR	Violet (D&C Violet n.2 - C.I. 60725) or undyed
SIZES AND STANDARD NEEDLES	EP from 0.2 to 7; USP from 10/0 to 5; Glicofil Lac is attached to AISI 300 or AISI 400 stainless steel atraumatic needles; "easy-detaching" needle type is also available to allow the detachment of the needle by a simple movement of the needle-holder
TENSILE STRENGTH	After 7 days from implant - 90% of the tensile strength remains, after 14 days - 75%, after 21 days - 50%, after 28 days - 25%.
ABSORPTION	Complete within 56-70 days.
INDICATIONS	In general surgery, thoracic surgery, orthopaedics, gynaecology, neurosurgery, in soft tissues approximation, when it is necessary an absorbable surgical suture.

THREAD LENGTH	SUFFIX
cm 45	AH
cm 60	AM
cm 70	AP
cm 75	AQ
cm 90	AT
cm 100	AV



Exclusive Distributor:

PT. LIMA MITRA ANDALAN
Ruko Grand Caman
Jl. Caman Raya No.54, Jati Bening, Bekasi
Indonesia 17412
Phone : +62 21 849 777 40 Fax: +62 21 849 0 1984
E-mail : info@limamitra.com web : www.limamitra.com

GLICOFIL® LAC

ABSORBABLE BRAIDED



BENANG

PGLA - Polyglycolic Lactid Acid

Efektifitas

- ⊙ Dapat memegang jaringan selama masa penyembuhan (masa pemulihan jaringan 35 hari)
- ⊙ Simpul kuat, sehingga jaringan tidak mudah lepas Mudah digunakan
- ⊙ Lembut, karena terbuat dari struktur benang yang halus dengan pelapis yang tipis (kurang dari 0.5%), sehingga mudah saat menarik benang pada waktu penjahitan
- ⊙ Dengan teknologi laser, penyatuan benang dan jarum menjadi lebih sempurna, sehingga kuat dan tidak menimbulkan trauma jaringan
- ⊙ Akses langsung, hanya dengan 2 langkah Nyaman bagi pasien
- ⊙ Hasil penyembuhan luka operasi sangat baik
- ⊙ Benang tidak perlu diangkat, sehingga tidak sakit dan tidak membuat cemas pasien

Aman

- ⊙ Respon inflamasi minimal, karena penyerapan secara hidrolisis
- ⊙ Masa absorpsi hari ke 56-70

JARUM

Berkualitas tinggi

- ⊙ Terbuat dari stainless steel dengan kualitas AISI 400. Kualitas stainless steel setara dengan jarum benang bedah cardiovascular Jarum tidak mudah bengkok

Ketajaman

- ⊙ Tetap tajam walaupun dilakukan sampai dengan 10x penetrasi
- ⊙ Dilapisi silikon

Indikasi

- ⊙ Bedah Umum
- ⊙ Bedah Thorax
- ⊙ Bedah Digestive
- ⊙ Bedah Syaraf
- ⊙ Bedah Gynecology
- ⊙ Bedah Urology
- ⊙ Bedah Orthopedic

GLICOFIL LAC

Polyglactin 910/PGLA : Multifilament, Synthetic, Mid term Absorbable suture



AGO -NEEDLE-AGUJA mm.	INDICATION	COLOR	THREAD X CM	USP	CODE	PACKS	UNIT PRICE (PCS)
TAPER 1/2 CIRCLE							
 HR 12,5	Dura	Violet	70	6 - 0	FY355L	36	Contact us
	Urologi	Violet	70	5 - 0	FX365L	36	Contact us
 HR 17,4	Urologi	Violet	70	5 - 0	FX374L	36	Contact us
	Usus, Dura	Violet	70	4 - 0	FW375L	36	Contact us
 HR 20	Urologi	Violet	70	5 - 0	FX384L	36	Contact us
	Usus	Violet	70	4 - 0	FW381L	36	Contact us
		Violet	70	3 - 0	FV382L	36	Contact us
 HR 25,9	Usus	Violet	100	3 - 0	FV386LAV	36	Contact us
		Violet	100	2 - 0	FU395LAV	36	Contact us
 HR 31,7	Usus	Violet	75	3 - 0	FV802L	36	Contact us
		Violet	75	2 - 0	FU799L	36	Contact us
 HR 31,1	Usus	Violet	75	3 - 0	FV245L	36	Contact us
		Violet	75	2 - 0	FU246L	36	Contact us
 HR 36,4	Subcutis, Otot	Violet	100	2 - 0	FU440LAV	36	Contact us
		Violet	100	0	FS464LAV	36	Contact us
 HRG 39,9	Subcutis, Otot	Violet	100	2 - 0	FU800LAV	36	Contact us
		Violet	100	0	FS445LAV	36	Contact us
	Uterus, Fascia	Violet	100	1	FR458LAV	36	Contact us
 HRG 48	Uterus, Fascia	Violet	100	1	FR479LAV	36	Contact us
	Uterus	Violet	100	2	FQ478LAV	36	Contact us

CUTTING 3/8 CIRCLE							
 TS 13,1	Subcutis	Violet	45	6 - 0	FY485L	36	Contact us
		Violet	45	5 - 0	FX503L	36	Contact us
 TS 14,7	Kulit	Violet	45	4 - 0	FW500L	36	Contact us
	Subcutis	Violet	45	3 - 0	FV528L	36	Contact us
 TS 18,7	Kulit	Violet	70	4 - 0	FW515L	36	Contact us
	Subcutis	Violet	70	3 - 0	FV516L	36	Contact us
 TS 24,3	Subcutis	Violet	70	3 - 0	FV533L	36	Contact us
	Kulit	Violet	70	2 - 0	FU534L	36	Contact us
 TFP 71,2	Kulit	Undyed	75	6 - 0	FY484LWAQ	36	Contact us
	Subcutis	Undyed	75	5 - 0	FX487LWAQ	36	Contact us
 TFP 11,9	Kulit	Undyed	75	5 - 0	FX489LWAQ	36	Contact us
 TFP 11,3	Subcutis	Undyed	75	5 - 0	FX511LWAQ	36	Contact us
	Subcutis	Undyed	75	4 - 0	FW870LWAQ	36	Contact us
 TFP 16,7	Subcutis	Undyed	75	4 - 0	FW812LW	36	Contact us
	Subcutis	Undyed	75	3 - 0	FV526LWAQ	36	Contact us
 TFP 14,1	Kulit	Undyed	75	3 - 0	FV523LWAQ	36	Contact us

ROUND BLUNT 1/2 CIRCLE							
 HT 40	Hepar	Violet	100	1	FR830L AV	36	Contact us

TAPER 5/8 CIRCLE							
 FRG 26,7	Laparoscopy	Violet	70	2-0	FU429	36	Contact us

Harga berlaku mulai bulan Januari 2016

Harga belum termasuk PPN

Harga dapat berubah sewaktu - waktu tanpa pemberitahuan terlebih dahulu